



Daniel Endoscopic Forehead Lift Surgery Set

Category: Plastic Surgery Set

Model / SKU: SS-3016

Product Page: <https://hamsansurgical.com/product/daniel-endoscopic-forehead-lift-surgery-set/>



Product Overview

The Daniel Endoscopic Forehead Lift Surgery Set is a sixteen piece instrument set for endoscopic brow and forehead lift procedures. Comprising curved and straight elevators, an orbital rim dissector, nerve hook, grasping forceps and dissection scissors, the set supports subperiosteal and subgaleal flap elevation through minimal scalp incisions.

Constructed from ISO Grade Stainless Steel and EU Medical Device compliant, the set is fully reusable and autoclavable, arriving in a stainless steel sterilizing box for immediate reprocessing readiness. Complete set procurement from a single supplier reduces sourcing complexity for plastic surgery departments and distributors.

Built for the Procedure.

Subgaleal forehead dissection requires elevators with defined curvature at the orbital rim and temporal planes. Instrument geometry designed for this anatomical access shapes the set's paired curved elevators.

Established reprocessing standards require autoclave-durable construction for instruments reused across endoscopic brow lift cases. ISO Grade Stainless Steel maintains structural integrity through repeated sterilisation cycles in this set.

Technical Specifications

Specification	Details
Product Name	Daniel Endoscopic Forehead Lift Surgery Set
Model / SKU	SS-3016
Material	ISO Grade Stainless Steel
Instrument Count	16 pieces
Reusable / Single-use	Reusable
Sterilisation Method	Autoclavable
Case / Tray	Stainless steel sterilizing box included

Clinical Applications

- Endoscopic forehead lift: subperiosteal and subgaleal flap dissection and elevation
- Brow ptosis correction: controlled repositioning of the brow-frontalis complex
- Lateral orbital rim release: periosteal release using the dedicated orbital rim dissector
- Glabellar region dissection: nerve-protected access near the supraorbital and supratrochlear nerves

Key Features

- Complete 16-piece set: covers full endoscopic forehead lift instrumentation in one case
- Left and right curved elevators: for lateral temporal and orbital rim access
- Orbital rim dissector: dedicated periosteal release along the lateral orbital rim
- Suction elevator: combines tissue elevation with continuous field clearing
- Nerve hook and nerve dissector: dedicated tools for the supraorbital nerve region
- Curved grasping forceps, 6": left and right curved pattern for tissue and flap handling
- Straight and curved dissection scissors: for controlled subgaleal dissection
- Stainless steel sterilizing box: included for organized storage and reprocessing

Usage and Care Guidelines

1. Inspect each instrument before use: check for visible damage, joint function, and tip integrity.
2. Clean immediately after use: rinse under running water to remove biological material.
3. Decontaminate using a validated washer-disinfector or manual cleaning protocol.
4. Sterilise by steam autoclave, per facility protocol.
5. Lubricate hinged joints with instrument-grade lubricant after each sterilisation cycle.
6. Store in the stainless steel sterilizing box or a clean, dry, covered tray.
7. Do not use if any component shows signs of corrosion, deformation, or functional impairment.
8. Follow facility reprocessing SOP and applicable national or international standards.

Frequently Asked Questions

Q: What material is the Daniel Endoscopic Forehead Lift Surgery Set made from?

A: The set is constructed from ISO Grade Stainless Steel.

Q: Is the Daniel Endoscopic Forehead Lift Surgery Set autoclave-compatible?

A: Yes, all sixteen instruments are fully reusable and autoclavable.

Q: What procedures is the Daniel Endoscopic Forehead Lift Surgery Set used for?

A: The set is used in endoscopic forehead and brow lift procedures, including subperiosteal and subgaleal flap dissection.

HAMSAN SURGICAL

Precision Crafted. Surgeon Trusted.

Shop Now: <https://hamsansurgical.com/product/daniel-endoscopic-forehead-lift-surgery-set/>

Browse All Products: <https://hamsansurgical.com>

© HAMSAN SURGICAL 2026. All rights reserved. This document is for informational purposes only. Specifications are subject to change without prior notice.

Scan For Product Details

